

A-dec Inspire[®] 391 Treatment Console



Monitor and sub base cover shown are not included

A-dec Inspire Treatment Console Features:

- Telescoping monitor mount for optimal monitor positioning throughout treatment procedure
- Dual floating shelves allow versatile positioning for procedural tubs and ancillary devices
- Four-drawer module to maximize in-treatment room organization and storage
- Single drawer for spare instruments, combined with a USB-powered vented storage area for the treatment room computer
- Side-opening access to a single swing-out shelf, ideally positioned for the dental assistant
- Sub base ready for cabinet mounted unit (541pro/341pro/545sta)
- Tip-out organizer bins
- Quiet soft-close drawers and doors
- Task lighting matches color temperature of the A-dec LED dental light
- USA and Canada UL-approved wiring with hospital-grade GFCI tamper-resistant outlets
- Single duplex outlet and USB 2.0
- Glove dispenser
- Clock timer
- 8 glass shelves
- Choice of pass-through or non-pass-through upper storage

391 quartz option:

- 2 cm quartz countertop
- 30" (762 mm) quartz countertop height
- 72" (1,829 mm) cabinet height (does not include clearance for telescoping monitor mount)
- 78" (1,981 mm) total height including clearance for telescoping monitor mount
- 22" (559 mm) depth of cabinet
- 42.5" (1,080 mm) width of cabinet (includes edge of the countertop to edge of side-opening handle)

391 solid surface option:

- 1 cm solid surface countertop
- 29.37" (746 mm) solid surface countertop height
- 72" (1,829 mm) cabinet height (does not include clearance for telescoping monitor mount)
- 78" (1,981 mm) total height including clearance for telescoping monitor mount
- 22" (559 mm) depth of cabinet
- 42.5" (1,080 mm) width of cabinet (includes edge of the countertop to edge of side-opening handle)

To learn more about the A-dec Inspire 391 Treatment Console, visit a-dec.com/dental-cabinets/a-dec-inspire or contact your local authorized dealer.

85.1042.00/BB/3.5M/4-23/Rev A